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
```

CDS      join(5883..6237,10472..10632,14063..14180,16680..16757,
        16897..17044,26646..26799,27717..27875)
        /gene="TGFB1"
        /gene_synonym="CED; DPD1; IBDIMDE; LAP; TGF-beta1; TGFB;
        TGFbeta"
    
```

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NCBI Reference Sequences (RefSeq)

RefSeqs maintained independently of Annotated Genomes

These reference sequences exist independently of genome builds. [Explain](#)

Genomic

1. **NG_013364.1 RefSeqGene**

Range	5005..28604
Download	GenBank , FASTA , Sequence Viewer (Graphics)

mRNA and Protein(s)

1. [NM_000660.7](#) → [NP_000651.3](#) transforming growth factor beta-1 proprotein preproprotein
[See identical proteins and their annotated locations for NP_000651.3](#)

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```
CDS      879..2051
         /gene="TGFB1"
         /gene_synonym="CED; DPD1; IBDIMDE; LAP; TGF-beta1; TGFB;
         TGFBeta"
```

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NCBI Resources How To Sign in to NCBI

ORFfinder PubMed Search

COVID-19 is an emerging, rapidly evolving situation.
 Get the latest public health information from CDC: <https://www.coronavirus.gov>.
 Get the latest research from NIH: <https://www.nih.gov/coronavirus>.
 Find NCBI SARS-CoV-2 literature, sequence, and clinical content: <https://www.ncbi.nlm.nih.gov/sars-cov-2/>.

Open Reading Frame Finder

ORF finder searches for open reading frames (ORFs) in the DNA sequence you enter. The program returns the range of each ORF, along with its protein translation. Use ORF finder to search newly sequenced DNA for potential protein encoding segments, verify predicted protein using newly developed SMART BLAST or regular BLASTP.

This web version of the ORF finder is limited to the subrange of the query sequence up to 50 kb long. Stand-alone version, which doesn't have query sequence length limitation, is available for [Linux x64](#).

Examples (click to set values, then click Submit button):

- NC_011604 Salmonella enterica plasmid pWES-1; genetic code: 11; 'ATG' and alternative initiation codons; minimal ORF length: 300 nt
- NM_000059; start codon: 'ATG only'; minimal ORF length: 150 nt

Enter Query Sequence

Enter accession number, gi, or nucleotide sequence in FASTA format:

```
>hTGFb1
at gccgcctccc gggctgccc
  901 tgcctccgct gctgctaccg ctgctgtggc tactggtgct gacgcttggc cggccggccg
  961 cgggactatc cactgcaag actatcgaca tggagctggg gaagcgggag cgcctcgagg
 1021 ccatccggcg ccagatcctg tccaagctgc ggctgcccag cccccgagc cagggggagg
 1081 tggccggcgg cccgctgccc gaggccgtgc tcgcccgtga caacagcacc cgcgaccggg
 1141 tggccgggga gaggcagaaa cgggagcccg agcctgaggc cgactactac gccaaaggag
 1201 tcaccggcgt gtaaatgggt gaaaccacac agaaatcta tgacaagttc aagcagagta
 1261 cacacagcat atatatgttc ttcaacacat cagagctccg agaagcggta cctgaaccgg
 1321 tgtgtctctc ccgggcagag ctgctctctc tgaggctcaa gttaaaagtg gacgagcag
```

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Open Reading Frame Viewer Help

Sequence

ORFs found: 7 Genetic code: 1 Start codon: 'ATG' only

ORF1 0 100 150 200 250 300 350 400 450 500 550 600 650 700 750 800 850 900 950 1 K 1,050 1,100 1,173

ORFfinder_9.25.234138898

ORF7 ORF2 ORF1 ORF6 ORF3 ORF4 ORF5

1: 1..1.2K (1,173 nt) Tracks shown: 2/5

Six-frame translation...

ORF1 (390 aa) Display ORF as... Unmark

```
>|c1|ORF1
MPPSGRLRLLPLLLPLLLVLTGPRPAAGLSTCKTIDMELVKKRRIEAIK
GQLLSKLRLASPPSQEVPVPGPLPEAVLALYNSTRDVRVAGESAEPEPEPE
ADYYAKVTRVLNVEITHIELYDKFKQSTHSYIFNITSELRLEAVPEPVL
SRAELRLLRLLKLVQHVLYQKYSNHSRHYLSNRLAPSDSEPHLSFDV
TGVNRQMLSRGGEIEGFRLSAHCSDSRDNTLQVDINGFTTGRGRLATZ
HGMNRPFLLLMATPLERAQLQSSRRRALDNYCFSSTKEMCCVQRQLYZ
DFRXDLGHWKIHPEPKYHANFCLGCPYINSLDTQYSKVLALYNQHNPQA
SAAPCCVQALEPLTVIYVYVGRKPKVEQLSNMIVRSCCKC
```

Label	Strand	Frame	Start	Stop	Length (nt aa)
ORF1	+	1	<1	1173	1173 390
ORF4	-	1	1038	586	453 150
ORF6	-	2	749	360	390 129
ORF3	+	2	953	1147	195 64

Mark subset... Marked: 1 Download marked set as CDS FASTA

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The screenshot shows the 'PrimerQuest Tool' interface. At the top, there are buttons for 'ASSAY DESIGN' and 'RESULTS'. On the right, there are buttons for 'HELP' and 'ABOUT'. The main section is titled 'Sequence Entry' and has a dropdown menu set to 'Enter sequence(s) manually'. Below this is a text box containing the DNA sequence from the previous image. To the right of the text box is a 'Hints' section with the following instructions:

- Paste your sequence into the Textbox
- Add up to 50 sequences in FASTA format
- Sequence length must be greater than 80 bases

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<input type="checkbox"/> Set 1 hTGFb1						
Amplicon Length: 1173 View Assay Details	Forward	Start	Stop	Length	Tm	GC%
		1	18	17	71	76.5
	Reverse	1153	1174	21	69	61.9
<input type="checkbox"/> Set 2 hTGFb1						
Amplicon Length: 1173 View Assay Details	Forward	1	18	17	71	76.5
	Reverse	1152	1174	22	71	63.6
<input type="checkbox"/> Set 3 hTGFb1						
Amplicon Length: 1173 View Assay Details	Forward	1	19	18	72	77.8
	Reverse	1153	1174	21	69	61.9
<input type="checkbox"/> Set 4 hTGFb1						
Amplicon Length: 1173 View Assay Details	Forward	1	19	18	72	77.8
	Reverse	1152	1174	22	71	63.6
<input type="checkbox"/> Set 5 hTGFb1						
Amplicon Length: 1173 View Assay Details	Forward	1	18	17	71	76.5
	Reverse	1154	1174	20	67	60

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Vary Reverse Primer Location
Off
Low
Med
High

Set Design Parameters for...

- General PCR (Primers Only)
- qPCR (2 Primers + Probe)
- qPCR Intercalating Dyes (Primers only)
- Forward Only - Sequencing
- Reverse Only - Sequencing

PrimerQuest Tool

ASSAY DESIGN RESULTS HELP ABOUT

Assay Set Locations for hTGFb1

ADD TO ORDER DOWNLOAD ASSAYS CUSTOMIZE ASSAY DESIGN SHOW CUSTOM TARGET REGIONS

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